# Finnpipette®

# Fixed Volume


Instructions for Use  
Bedienungsanleitung  
Guide d'utilisation et d'entretien  
Instrucciones de uso


## Thermo

ELECTRON CORPORATION


Thermo Electron Oy  
P.O.Box 100, Fin-01621 Vantaa, Finland  
Tel. +358-9-329 100, fax -358-9-3291 0414  
[www.thermo.com/finnpipette](http://www.thermo.com/finnpipette)


1


3


5


4


6


2


8


10


9


11

7

12

13

## CONTENTS ENGLISH

PRODUCT DESCRIPTION	4
PIPETTE OPERATION	4
PIPETTING TECHNIQUES	5
CALIBRATION	6
MAINTENANCE	8
TROUBLE SHOOTING	10
PACKAGE	10
SPARE PARTS	38

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## PRODUCT DESCRIPTION

The Finnpiquette Fixed Volume is an autoclavable digital pipette. It operates on the air displacement principle (i.e. an air interface) and uses detachable, disposable tips. The delivery volume is displayed on a readout window on the handle.

The sixteen different models of Finnpiquette Fixed Volume pipettes cover a volume range from 1 µl to 10 ml.

Order No.	Volume Range	Finntip
4501000	1 µl	10
4501010	2 µl	10
4501020	5 µl	10
4501030	10 µl	10
4501040	20 µl	250, 200 Ext
4501050	25 µl	250, 200 Ext
4501060	50 µl	250, 300, 200 Ext
4501070	100 µl	250, 300, 200 Ext
4501080	200 µl	250, 300, 200 Ext
4501090	250 µl	1000
4501100	500 µl	1000
4501110	1000 µl	1000
4501120	2 ml	5 ml
4501130	3 ml	5 ml
4501140	5 ml	5 ml
4501150	10 ml	10 ml

### 1 DIGITAL DISPLAY

The delivery volume is clearly indicated in the large digital display on the handle.

### RAW MATERIALS

The Finnpiquette Fixed Volume is made of mechanically durable and chemically resistant materials which allow repeated autoclaving of the complete pipette at 121°C.

### DESCRIPTION OF TIPS

Finntips are recommended for use with the Finnpiquette Fixed Volume.

They are made of natural colour polypropylene, generally regarded as the only contamination free material suitable for tips. Finntips are also autoclavable (121°C).

## PIPETTE OPERATION

### 2 TIP EJECTION

To help eliminate the risk of contamination, each pipette is fitted with a tip ejector system. The tip ejector system consists of a soft-touch tip ejector and specially designed gearing mechanism. To release the tip, point the pipette at suitable waste receptacle and press the tip ejector with your thumb.

**SAFETY LABEL**

12

You can mark the pipette application your initials the calibration date, etc. on the safety label. Remove the clear plastic window on the edge closest to the push button (use the service tool that comes with the pipette, or a screwdriver). Mark the adhesive label with a felt-tipped or other pen and snap the window back in place.

**SHELF HANGER**

13

You can attach the pipette shelf hanger on a counter, pipette stand or anywhere where you want to hang your pipette.

Clean the area where you plan to attach the shelf hanger. Apply two stickers to the underside of the shelf hanger. Press the shelf hanger firmly into place -- on a shelf, countertop or pipette stand. To use, hang the grippy finger rest on the shelf hanger.

**PIPETTING TECHNIQUES**

Push and release the push button slowly at all times, particularly when working with high viscosity liquids. Never allow the push button to snap back. Make sure that the tip is firmly attached to the tip cone. Check for foreign particles in the tip.

Before you begin your actual pipetting work, fill and empty the tip 2-3 times with the solution that you will be pipetting. Hold the pipette in an upright position while aspirating liquid. The grippy should rest on your index finger. Make sure that the tips, pipette and solution are at the same temperature.

**Figures 3-6:****A = Ready position****B = First stop****C = Second stop****FORWARD TECHNIQUE**

3

Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the push button to the first stop.
2. Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm and slowly release the push button. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
3. Deliver the liquid by gently depressing the push button to the first stop. After a delay of about one second, continue to depress the push button all the way to the second stop. This action will empty the tip.
4. Release the push button to the ready position.

If necessary, change the tip and continue pipetting.

**REVERSE TECHNIQUE**

4

The reverse technique is suitable for dispensing liquids that have a high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes.

Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the push button all the way to the second stop.
2. Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm, and slowly release the push button. This action will fill the tip. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.

3. Deliver the preset volume by gently depressing the push button to the first stop. Hold the push button at the first stop. Some liquid will remain in the tip and this should not be included in the delivery.
4. The remaining liquid should either be discarded with the tip or pipetted back into the container.

## 5 REPETITIVE TECHNIQUE

The repetitive technique offers a rapid and simple procedure for repeated delivery of the same volume. Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the push button all the way to the second stop.
2. Dip the tip under the surface of the liquid in the reservoir to a depth of about 1 cm, and slowly release the push button. This action will fill the tip. Withdraw the tip from the liquid touching against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the push button to the first stop. Hold the push button at the first stop. Some liquid will remain in the tip and this should not be included in the delivery.
4. Continue pipetting by repeating steps 2 and 3.

## 6 PIPETTING WHOLE BLOOD

(deproteinization in blood glucose determination, for example)

Use steps 1 and 2 of the forward technique to fill the tip with blood.

Wipe the tip carefully with a dry clean tissue.

1. Immerse the tip into the reagent and depress the push button to the first stop, making sure the tip is well below the surface.
2. Release the push button slowly to the ready position. This will fill the tip. Keep the tip in the solution.
3. Depress the push button to the first stop and release slowly. Keep repeating this procedure until the interior wall of the tip is clear.
4. Finally, depress the push button all the way to the second stop to completely empty the tip.

## CALIBRATION

All Finnpiettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. The pipettes are constructed to permit re-adjustment for liquids of different temperature and viscosity.

### DEVICE REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range	readable graduation
--------------	---------------------

under 10 $\mu$ l	0.00 1 mg
------------------	-----------

10-100 $\mu$ l	0.01 mg
----------------	---------

above 100 $\mu$ l	0.1 mg
-------------------	--------

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ( $\pm 0.5^{\circ}\text{C}$ ) temperature of water, pipette and air between  $15^{\circ}\text{C}$  to  $30^{\circ}\text{C}$ . The relative humidity must be above 50%. Especially with

volumes under 50 µl the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

## PROCEDURE TO CHECK CALIBRATION

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume. A new tip is first pre-wetted 3-5 times and a series of ten pipettings is done with both volumes. A pipette is always adjusted for delivery (Ex) of the selected volume.

### Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the accuracy (A) and precision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the calculated results are within the selected limits, the adjustment of the pipette is correct.

TABLE 1: Maximum permissible errors according ISO8655

Volume µl	Accuracy		Precision	
	µl	%	s.d.* µl	cv%
1	±0.04	±4.0	0.04	4.0
2	±0.05	±2.5	0.04	2.0
5	±0.07	±1.4	0.07	1.4
10	±0.09	±0.9	0.08	0.8
20	±0.12	±0.6	0.10	0.5
25	±0.15	±0.6	0.13	0.5
50	±0.3	±0.6	0.2	0.4
100	±0.4	±0.4	0.3	0.3
200	±0.8	±0.4	0.6	0.3
250	±1.0	±0.4	0.8	0.3
500	±1.5	±0.3	1.5	0.3
1000	±3.0	±0.3	0.3	0.3
2000	±6.0	±0.3	4.0	0.2
3000	±9.0	±0.3	6.0	0.2
5000	±15.0	±0.3	10.0	0.2
10 ml	±30.0	±0.3	20.0	0.2

## ADJUSTMENT

Adjustment is done with the service tool.

1. Place the service tool into the openings of the calibration nut at the top of the handle.
2. Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.
3. After adjustment check the calibration according to the instructions above.

## FORMULAS FOR CALCULATING RESULTS

### Conversion of mass to volume

$$V = (w + e) \times Z$$

V = volume (µl)

w = weight (mg)

e = evaporation loss (mg)

Z = conversion factor for mg/µl conversion

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6 mg = 0.2 mg/s).

Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2 mg (10 s x 0.2 mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary.

The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032 µl/mg at 22°C and 95 kPa. See the conversion table on page 36.

**Accuracy (systematic error)**

Accuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0$$

$\frac{A}{\bar{V}}$	= accuracy
$\bar{V}$	= mean volume
$V_0$	= nominal volume

Accuracy can be expressed as a relative value: **A% = 100% x A / V<sub>0</sub>**

**Precision (random error)**

Precision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$S = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

s	= standards deviation
$\bar{v}$	= mean volume
n	= number of measurements

Standard deviation can be expressed as a relative value (CV) **CV = 100% x S /  $\bar{V}$**

**MAINTENANCE**

When the Finnpiquette Fixed Volume is not in use, make sure it is stored in an upright position. We recommend a Finnpiquette stand for this purpose.

**SHORT-TERM CHECKING**

The pipette should be checked at the beginning of each day for dust and dirt on the outside surfaces of the pipette.

Particular attention should be paid to the tip cone. No other solvents except 70 % ethanol should be used to clean the pipette.

**LONG-TERM MAINTENANCE**

If the pipette is used daily it should be checked every three months. The servicing procedure starts with the disassembly of the pipette.

**DISASSEMBLING 1-40 µl PIPETTES**

- 8 1. Press the tip ejector.
2. Insert the maintenance pliers under the ejector bar to release the tip ejector.
- 9 3. Remove the tip cone by pressing with maintenance pliers.
4. Pull out the piston and the spring.
5. Keep the tip cone vertically and push out with piston the rest of the piston assembly. Then keep the tip cone upside down and tap all parts from tip cone. Remember keep all parts in order on table for reassembly.
6. Clean the piston, the piston spring and the O-rings with a dry napless cloth.
7. Check the tip cone for foreign particles.
8. Grease the cleaned parts with the lubricant that comes with the pipette.
9. Reassemble the pipette components.

**All 1-40 µl:** First, slide the spring 14, spring support 15 and tube 16 back on the piston. Compress the spring with fingers by pressing piston and spring support 15 against each other.

**11-40 µl:** Slide bigger O-ring 17, smaller O-ring 18, spring support 19 (sharp edges against spring) and small spring 20 on the piston.

**3-10 µl:** First slide O-ring tube 17 (larger hole first), bigger O-ring 18, smaller O-ring 19 and O-ring support 20 on the piston. Then slide small spring 21, spring support 22 (sharp edges against spring) and O-ring 22 on the O-ring support 20.

**1-2 µl:** First slide O-ring tube 17 (larger hole first) and sealing combination 18 on the piston. Then slide small spring 19, spring support 20 (sharp edges against spring) and O-ring 21 on the sealing combination 18.

**All 1-40 µl:** Carefully slide the entire assembly into the tip cone and release your fingers. **14**

10. With the push button depressed all the way carefully attach the tip cone to the handle so that the adapter opening is on the tip ejector side. Do not bend the thin piston wire when assembling. Press in the snap joints.
11. Assemble the tip ejector and check the calibration according to the instructions.

### DISASSEMBLING 41-1000 µl PIPETTES

1. Press the tip ejector. **8**
2. Insert the pliers under the ejector bar to release the tip ejector.
3. Remove the tip cone using the maintenance pliers. **9**
4. Pull out the piston.
5. Remove the O-ring, O-ring support and spring from the tip cone.
6. Clean the piston, the piston spring and the O-ring with a dry napless cloth.
7. Check the cylinder for foreign particles.
8. Grease the cleaned parts with the lubricant that comes with the pipette.
9. Slide parts over the piston pressing down the large spring. Attach the tip cone to the handle so that the adapter opening is on the tip ejector side, and press in the snap joints.
10. Check the calibration according to the instructions.

### DISASSEMBLING 2-10 ml PIPETTES

1. Press the tip ejector. **8**
2. Insert the pliers under the ejector bar to release the tip ejector.
3. Remove the part 2 from part 1 of the tip ejector using the maintenance pliers to release the snap joint. **10**
4. Remove the cylinder by pressing part 1 of the tip ejector firmly towards the cylinder. This action releases the snap joint so you can remove the cylinder. **11**
5. Clean the O-ring and cylinder. Regrease the O-ring and the cylinder.
6. Assemble the parts in the opposite order of disassembly. All joints are snap fit and can be pushed together by hand. Be careful not to bend the pipette during assembly because this could damage the snap joints or the piston.
7. Check the calibration according to the instructions.

### STERILIZATION

The entire pipette can be sterilized by autoclaving it at 121°C (252°F) (minimum 20 minutes). No special preparations are needed for autoclaving. You can use steam sterilization bags if needed.

After autoclaving the pipette must be cooled to room temperature for at least two hours. Before pipetting, make sure that the pipette is dry. We recommend that you check the calibration of 1-1000 µl pipettes after every 25th sterilization cycle and of 2-10 ml pipettes after every 10th sterilization cycle.

**CAUTION!**

*The Finnpiquette is designed to allow easy in-lab service. If you would prefer to have us or your local representative service your pipette, please make sure that the pipette has been decontaminated before you send it to us.*

*Please note that the postal authorities in your country may prohibit or restrict the shipment of contaminated material by mail.*

**TROUBLE SHOOTING**

The table below lists possible problems and their solutions.

Defect	Possible reason	Solution
Leakage	Tip incorrectly attached	Attach firmly
	Foreign particles between tip and tip cone	Clean tip cones attach new tips
	Foreign particles between the piston, the O-ring and the cylinder	Clean and grease O-ring and cylinder.
	Insufficient amount of grease on cylinder and O-ring	Grease accordingly
Inaccurate dispensing	O-ring damaged	Change the O-ring
	Incorrect operation	Follow instructions carefully
	Tip incorrectly attached	Attach firmly
Inaccurate dispensing with certain liquids	Calibration altered: caused by misuse, for example	Recalibrate according to instructions
	Unsuitable calibration. High viscosity liquids may require recalibration.	Recalibrate with the liquids in question.

**PACKAGE**

The Finnpiquette Fixed Volume is shipped in a specially designed package containing the following items:

- |                       |                       |                            |
|-----------------------|-----------------------|----------------------------|
| 1. The Finnpiquette   | 4. Finntip sample     | 7. Calibration certificate |
| 2. Service tool       | 5. Tube of grease     | 8. Shelf hanger            |
| 3. Maintenance pliers | 6. Instruction manual | 9. Two stickers            |

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Product specifications are subject to change without prior notice. Finnpiquette® and Finntip® are registered trademarks of Thermo Electron Oy.