



pipetman **NEO**®

The Pipetting Comfort with Improved Comfort!



EN USER'S GUIDE



IMPORTANT NOTICE!

*Do not lubricate this pipette.
The use of lubricant cancels
the warranty of this pipette.*

<i>Contents</i>	<i>Page</i>
1- INTRODUCTION	2
2- PARTS CHECK LIST	3
3- DESCRIPTION	3
4- SETTING THE VOLUME	4
5- PIPETTING	5
6- GENERAL GUIDELINES FOR GOOD PIPETTING	9
7- ACCESSORIES	10
8- GLP FEATURES	11
9- TROUBLESHOOTING	11
10- LEAK TEST	12
11- MAINTENANCE	13
12- CLEANING AND DECONTAMINATION	17
13- SPECIFICATIONS	21
14- SPARE PARTS	23
EC DECLARATION OF CONFORMITY	26
WARRANTY	27

1 - INTRODUCTION

PIPETMAN Neo is a fully adjustable air displacement pipette and is used with disposable tips.

PIPETMAN Neo is an evolution of the PIPETMAN Classic to address the growing susceptibility of pipette users to repetitive strain injuries. The implementation of new features to improve pipetting comfort has no impact on the legendary robustness, accuracy and precision of PIPETMAN Classic.


- ▶ Six single channel models cover a volume range from 0.2 μ L to 1000 μ L.
- ▶ Two 8 and two 12 multichannel models cover a volume range from 2 μ L to 200 μ L.


2 - PARTS CHECK LIST

Just take a moment to verify that the following items are present:

- PIPETMAN Neo,
- User's Guide,
- Safety bag,
- Certificate of conformity (including bar-code sticker).

3 - DESCRIPTION

 Please refer to the following chapters for a full description of the different parts and functions of the pipette.





4 - SETTING THE VOLUME

The volume of liquid to be aspirated is set using the volumeter.

Model	Color of volumeter numbers		
	Black	Red	Increment
P2N	μL	0.01 μL	0.002 μL
P10N to P20N	μL	0.1 μL	0.02 μL
P100N-P200N	μL	-	0.2 μL
P1000N	0.01 mL	mL	0.002 mL

P2N 1.25 μL	P10N 7.5 μL	P20N 12.5 μL
P100N 75 μL	P200N 125 μL	P1000N 0.75 mL

The dials are colored either black or red to indicate the position of the decimal point, depending on the model (see examples).

The volume is set by turning the thumbwheel or the push-button. The push-button makes it easier and quicker to set volumes, especially when wearing gloves. The thumbwheel may be turned using only one hand to slowly reach the required setting.




To obtain maximum accuracy when setting the volume, proceed as follows:

- ▶ when **decreasing** the volume setting, slowly reach the required setting, making sure not to overshoot the mark.
- ▶ when **increasing** the volume setting, pass the required value by 1/3 of a turn and then slowly decrease to reach the volume, making sure not to overshoot the mark.

5 - PIPETTING


For optimum performance, use of PIPETMAN DIAMOND Tips with your PIPETMAN Neo is strongly recommended. These Tips, made from pure polypropylene have the Gilson logo engraved on their collar, ensuring that you have a genuine Gilson product. Plastic tips are for a single application – they must not be cleaned for reuse.

 *PIPETMAN Neo can also be used with the main tip brands.*

Fitting the tips

Single channel models:

To fit a new PIPETMAN DIAMOND Tip, push the tip-holder into the tip using a slight twisting motion to ensure a firm, airtight seal.

 *For the P2N and P10N models, a dual-position adapter (plastic) is required to fit DL10 tips (long tips) or D10 tips (short tips). The metallic rod of the tip-ejector is shaped so that the adapter may be clipped to it in either of two positions.*

P2N and P10N models are delivered with the adapter in place, positioned in the longer slot, ready to use DL10 tips. When D10 tips (which are shorter) are used, the adapter must be repositioned in the shorter slot as follows:

- 1 Pull the adapter down from the metallic rod.
- 2 Turn the adapter through 180°C.
- 3 Refit the adapter so that the end of the metallic rod engages the shorter slot of the adapter.



Multichannel models:

PIPETMAN DIAMOND Tip compatibility

8x20 DL10, DFL10, DF30, D200
12x20

8x200 D200, D300, DF100, DF200, DF300
12x200



PIPETMAN DIAMOND Tip: TIPACK™ and TOWERPACK™ – ROCKY RACK™

PIPETMAN DIAMOND Tips are best fitted from the patented ROCKY RACK available only in our TIPACKS and TOWER PACKS.

ROCKY RACK is the dome-shaped part of the pack that contains the tips. ROCKY RACK makes it easy to securely fit the tips to a multichannel pipette, ensuring an airtight seal on all channels without the need to use undue pressure or to touch the tips.

The patented ROCKYRACK system available on TOWERPACK and TIPACK makes it extremely easy to fit tips on a multichannel pipette. Tips will not fall off nor will they have to be manually positioned.




Pre-rinse the tips

Some liquids (e.g. protein-containing solutions and organic solvents) can leave a film of liquid on the inside wall of the tip ; pre-rinse the tip to minimize any errors that may be related to this phenomenon.

Pre-rinsing consists of aspirating the first volume of liquid and then dispensing it back into the same vessel (or to waste). Subsequent volumes that you pipette will have levels of accuracy and precision within specifications.

Aspirate

- 1 Press the push-button to the **first stop** (this corresponds to the set volume of liquid).
- 2 Hold the pipette vertically and immerse the tip in the liquid (see immersion depth table, page 10). Release the push-button slowly and smoothly (to **top** position) to aspirate the set volume of liquid. Wait one second (time depends on model, see table) ; then withdraw the pipette-tip from the liquid. You may wipe any droplets away from the outside of the tip using a medical wipe, however if you do so take care to avoid touching the tip's orifice.



Dispense

- 1 Place the end of the tip against the inside wall of the recipient vessel (at an angle of 10° to 40°).
- 2 Press the push-button slowly and smoothly to the **first stop**.

Wait for at least a second, then press the push-button to the **second stop** to expel any residual liquid from the tip. Keep the push-button pressed fully down and (while removing the pipette) draw the tip along the inside surface of the vessel.

- 3 Release the push-button, smoothly. Eject the tip by pressing firmly on the tip-ejector button.

6-GENERAL GUIDELINES FOR GOOD PIPETTING

- 1 Make sure that you operate the push-button slowly and smoothly.
- 2 When aspirating, keep the tip at a constant depth below the surface of the liquid (refer to the table).
- 3 Change the tip before aspirating a different liquid, sample, or reagent.

Table - Immersion Depth and Wait Time

Model	Immersion Depth (mm)	Wait Time (seconds)
P2N	1	1
P10N	1	1
P20N	2-3	1
P100N	2-4	1
P200N	2-4	1
P1000N	2-4	2-3
8x20, 12x20	2-3	1
8x200, 12x200	2-3	1

- 4 Change the tip if a droplet remains at the end of the tip from the previous pipetting operation.
- 5 Each new tip should be pre-rinsed with the liquid to be pipetted.
- 6 Liquid should never enter the tip-holder; to prevent this:
 - press and release the push-button slowly and smoothly,
 - never turn the pipette upside down,
 - never lay the pipette on its side when there is liquid in the tip.
- 7 If you use the same tip with a higher volume, pre-rinse the tip.
- 8 For volatile solvents you should saturate the air-cushion of your pipette by aspirating and dispensing the solvent repeatedly before aspirating the sample.
- 9 When the temperature of the liquid is different from the ambient temperature, pre-rinse the tip several times before use.
- 10 You may remove the tip-ejector (see Chapter 11 - Maintenance) to aspirate from very narrow tubes.
- 11 After pipetting acids or other corrosive liquids that emit vapors, remove the tip-holder, rinse the

piston, O-ring and seal with distilled water. For the model P1000N, by using a specific tip holder equipped with a filter, you can increase the lifetime of the piston (see Chapter 7 - Accessories).

- Do not pipette liquids having temperatures above 70 °C or below 4 °C. The pipette can be used between + 4 °C and + 40 °C but the specifications may vary according to the temperature (refer to the ISO 8655-2 standard for conditions of use).

7 - ACCESSORIES

To make pipetting more comfortable and more secure, Gilson has developed several accessories:

- To avoid the possibility of liquid running back into the pipette, store the pipette vertically.

CARROUSEL™ Pipette stand (7 pipettes)	F161401
TRIO™ stand (3 pipettes)	F161405
SINGLE™ pipette holder	F161406

- To identify or personalize your pipette, COLORIS™ clips are available:

COLORIS™ clips (mixed colors set of 10)	F161301
COLORIS™ clips (red, set of 10)	F161302
COLORIS™ clips (yellow, set of 10)	F161303
COLORIS™ clips (green, set of 10)	F161304
COLORIS™ clips (blue, set of 10)	F161305
COLORIS™ clips (white, set of 10)	F161306

- With THE JIMMY™, handsfree microtube opener, you can open both snap-cap and screw-cap microtubes.

THE JIMMY™ (set of 3)	F144983
-----------------------	---------

- To protect the piston when pipetting corrosive liquids, you can use a specific tip holder and filter for the model P1000N:

Corrosion protection kit (tip holder + a bag of 10 filters)	F144570
---	---------

8 - GLP FEATURES


The **Serial Number** is engraved on the body of the pipette. It provides unique identification of your pipette and the date of manufacture; Ex : AA 10369

The **Bar Code** on the box and the certificate of conformity provide traceability of your pipette.



9 - TROUBLESHOOTING

A quick inspection of the pipette may help you to detect a problem.

 You may download from the Gilson website (www.gilson.com) the "2 minute inspection", which shows how to perform a quick diagnosis of your pipette.



Before returning any pipette to your local Gilson Service Center, ensure that it is completely free of chemical, biological, or radioactive contamination. Refer to Chapter 12 - Cleaning and Decontamination. Please use the included safety bag to return the pipette to your local Gilson Service Center.

The following table may help you to identify and correct the problem you might encounter.

Symptom	Possible Cause	Action
Pipette is leaking sample	Damaged tip-holder	Replace the tip-holder*
	Worn O-ring or seal	Replace both parts*
Pipette won't aspirate	Worn O-ring	Replace both parts*
	Damaged tip-holder	Replace the tip-holder*
	Connecting nut is loose	Tighten connecting nut
	Damaged or corroded piston	Return pipette to supplier
	Improper repair or assembly	See Chapter 11 - Maintenance

Symptom	Possible Cause	Action
Pipette is not precise	Tip-holder is loose	Tighten connecting nut*
	Connecting nut is loose	Tighten connecting nut
	Incorrect operator technique	Operator training
	Damaged or corroded piston(s)	Return pipette to supplier
	Damaged tip-holder(s)	Replace the tip-holder*
	Worn O-ring or seal	Replace both parts*
Tips fall off or do not fit correctly	Low quality tips	Use PIPETMAN DIAMOND tips
	Damaged tip-holder(s)	Replace the tip-holder*
	Damaged tip-ejector	Replace tip-ejector
	Ejector spacer is missing	Mount the spacer on the tip- ejector
	The ejector spacer is damaged	Replace the ejector spacer
	The tip-ejector is loose	Assemble the tip-ejector properly
	The ejector lock is misaligned	Align the ejector lock

* Not applicable for multichannel models

However, if you can't solve the problem, contact your Gilson representative.

10 - LEAK TEST

This test may be performed at any time to check that the pipette does not leak, especially after performing a maintenance or decontamination procedure. If a pipette fails this test, replace the O-ring and seal. After making sure that the pipette is correctly reassembled, repeat this test (not applicable to multichannel pipettes).

For the P2N to P200N, 8x20, 12x20, 8x200 and 12x200 models:

- 1 Fit a PIPETMAN DIAMOND Tip.
- 2 Set the pipette to the maximum volume given in the specifications, and pre-rinse.
- 3 Aspirate the set volume from a beaker of distilled water.
- 4 Maintain the pipette in the vertical position and wait for 20 seconds.
- 5 If a water droplet appears at the end of the tip there is a leak.
- 6 If you see no droplet, re-immers the tip below the surface of water.

- 7 The water level inside the tip should remain constant; if the level goes down there is a leak.

For the P1000N model:

- 1 Fit a PIPETMAN DIAMOND tip.
- 2 Set the pipette to the maximum volume given in the specifications.
- 3 Aspirate the set volume from a beaker of distilled water.
- 4 Maintain the pipette in the vertical position and wait for 20 seconds.
- 5 If a water droplet appears at the end of the tip, there is a leak.

11 - MAINTENANCE

Routine maintenance will help keep your pipette in good condition, ensuring a continued high level of performance. Maintenance is limited to cleaning or autoclaving the parts specified under Chapter 12 - Cleaning and Decontamination or to replacing the push-button, connecting nut, tip-ejector, tip-holder, seal and O-ring.

PIPETMAN Neo P2N and P10N should not be disassembled, so you may only replace the push-button, tip-ejector, dual position tip-ejector and its adapter. With these pipettes if the tip-holder is damaged, the piston may also be damaged.

For multichannel models, the lower part should not be disassembled: only the push-button, the connecting nut and the tip-ejector can be replaced.



After replacing any parts you should verify the performance of your pipette following the verification procedure available on the Gilson website (www.gilson.com). If the pipette needs to be readjusted, please contact your local Gilson authorized Service Center.

Single channel models

Changing the Tip-ejector, to remove

- 1 Push the ejection button.
- 2 Push laterally the tip-ejector.
- 3 Slide and remove the tip-ejector.



To refit

- 1 Push the ejection button.
- 2 Slide the tip-ejector along the tip-holder.
- 3 Clip the tip-ejector on the body of the pipette.


Changing the Tip-holder – no tools required

- 1 Remove the tip-ejector (see above).
- 2 Unscrew the connecting nut by turning it counter-clockwise.
- 3 Carefully separate the lower and upper parts.
- 4 Remove the piston assembly, O-ring and seal.
- 5 Clean, autoclave, or replace the tip-holder.
- 6 Reassemble the pipette (refer to the figure, p 15).
- 7 Tighten the connecting nut (turn clockwise).
- 8 Refit the tip-ejector (see above).

Servicing the Piston Assembly

You may remove the piston assembly for cleaning purposes only. If the piston assembly is changed, the pipette must be adjusted and calibrated in a Gilson authorized Service Center.

- 1 Remove the tip-ejector (see above).
- 2 Unscrew the connecting nut by turning it counter-clockwise.
- 3 Carefully separate the lower and upper parts.
- 4 Remove the piston assembly, O-ring and seal.
- 5 Clean and decontaminate the piston assembly.



The piston assembly must not be autoclaved.

- ⑥ Reassemble the pipette (refer to the figure, page 15).
- ⑦ Tighten the connecting nut (turn clockwise).
- ⑧ Refit the tip-ejector (see above).

Changing the O-ring

The O-ring and seal are on the piston; **they must not be autoclaved**, if worn or damaged in any way (chemical or mechanical), they must be replaced.

The dimensions of the O-ring vary depending on the model of pipette.

- ① Remove the tip-ejector (see above).
- ② Unscrew the connecting nut by turning it counter-clockwise.
- ③ Carefully separate the lower and upper parts.
- ④ Remove the piston assembly, O-ring and seal.
- ⑤ Clean or replace the seal then the O-ring.
- ⑥ Reassemble the pipette (refer to the figure, page 15).
- ⑦ Tighten the connecting nut (turn clockwise).
- ⑧ Refit the tip-ejector (see above).

Multichannel models

Changing the Tip-ejector

- ① To remove the tip-ejector, keep both ejector locks depressed.
- ② Pull the tip-ejector down.
- ③ To refit the tip-ejector, gently reinsert the tip-ejector vertically into the rails of the ejector support.
- ④ Pull lightly on the tip-ejector to check the position.



12 - CLEANING AND DECONTAMINATION

PIPETMAN Neo is designed so that the parts normally in contact with liquid contaminants, can easily be cleaned and decontaminated. However, because the models P2N and P10N contain miniaturized parts, it is best not to disassemble these pipettes yourself; please contact your local Gilson authorized Service Center.



*You may refer to the decontamination procedure available on the Gilson website (www.gilson.com). **Liquid must never enter the upper part (handle) of any pipette.***

Cleaning

The pipette must be cleaned, as described below, before it is decontaminated. Soap solution is recommended for cleaning PIPETMAN Neo.

Single channel models - External

- ① Remove the tip-ejector.
- ② Wipe the tip-ejector with a soft-cloth or lint-free tissue impregnated with soap solution.
- ③ Wipe the entire pipette with a soft-cloth or lint-free tissue impregnated with soap solution, to remove all dirty marks. If the pipette is very dirty, a brush with soft plastic bristles may be used.
- ④ Wipe the entire pipette and the tip-ejector with a soft cloth or lint-free tissue soaked with distilled water.
- ⑤ Refit the tip-ejector and allow the pipette to dry.

Single channel models - Internal

The following components **only** can be immersed in a cleaning solution: connecting nut, tip-ejector, tip-holder, piston assembly, seal and O-ring.

- ① Disassemble the pipette as described in the Chapter 11 - Maintenance.
- ② Set aside the upper part in a clean, dry place.
- ③ Clean the individual components of the lower part of the pipette using an ultrasonic bath (20 minutes at

50°C) or with a soft-cloth and brushes. Small round brushes with soft plastic bristles may be used to clean the interior of the tip-holder.

- ④ Rinse the individual components with distilled water.
- ⑤ Leave the parts to dry by evaporation or wipe them with a clean soft-cloth or lint-free tissue.
- ⑥ Reassemble the pipette as described in the Chapter 11 - Maintenance.

Multichannel models

The following components **only** can be immersed in a cleaning solution: tip-ejector, ejector locks and ejector spacer.

- ① Remove the tip-ejector and the ejector spacer.
- ② Immerse the tip-ejector, ejector locks and ejector spacer in the cleaning solution or wipe them with a soft-cloth or lint-free tissue impregnated with the cleaning solution.
- ③ Rinse the components with distilled water.
- ④ Wipe the entire pipette with a soft cloth or lint-free tissue impregnated with the cleaning solution.
- ⑤ Wipe it with distilled water.
- ⑥ Leave the parts to dry by evaporation or wipe them with a clean soft-cloth or lint-free tissue.
- ⑦ Refit the tip-ejector as described in page 16 (Changing the tip-ejector).

Autoclaving

Single channel models

The upper part (body) and the piston assembly of the pipette are **not** autoclavable. **Only** the following parts may be autoclaved: tip-ejector, tip-holder and connecting nut. The O-ring and seal are **not** autoclavable; they may be cleaned or replaced with the one specified in Chapter 14 - Spare Parts.

- ① Clean the parts to be autoclaved, especially the tip-holder.
- ② Put the parts in an autoclaving sack.
- ③ Autoclave for 20 minutes at 121°C, 0.1 MPa.

- ④ Check that the parts are dry before re-assembling the pipette.
- ⑤ Set the pipette aside to stabilize at room temperature.

Multichannel models

The complete lower part can be autoclaved (20 minutes at 121°C, 0.1 MPa).

The lower part has to be removed from the body.

To remove:

- ① Keep the tip-ejector button depressed,
- ② Turn the ejector clip counter-clockwise and separate it from the handle.
- ③ Unscrew the connecting nut turning counter-clockwise.
- ④ Press push-button to extract the head.

To autoclave:

- ① Clean the parts to be autoclaved, especially the tip-holders and the ejector spacer.
- ② Put the lower part in an autoclaving sack.
- ③ Autoclave for 20 minutes at 121°C, 0.1 MPa.
- ④ Check that the autoclaved part is dry before re-assembling the pipette.
- ⑤ Set the pipette aside to stabilize at room temperature.



To fit the lower part on the body,

- 1 Press on the push-button.
- 2 Insert the lower part into the ejection rod. Release the push-button.
- 3 Screw the connecting nut turning clockwise.
- 4 Press on the tip-ejector button.
- 5 Insert the ejector clip and release the tip-ejector button.



Chemical Decontamination

You may choose to decontaminate your pipette chemically, in accordance with your own procedures. Whatever decontaminant you use, check with the supplier of the decontaminant that it is compatible with stainless steel and the plastics used in the construction of the pipette: PA (Polyamide), PBT (Polybutylene Terephthalate), PC (Polycarbonate), PC/PBT (Polycarbonate/ Polybutylene Terephthalate), POM (Polyoxymethylene), or PVDF (Polyvinylidene Fluoride).

Single and Multichannel - Upper Part (handle)

- 1 Wipe the upper part (handle) of the pipette with a soft-cloth or lint-free tissue impregnated with the chosen decontaminant.
- 2 Wipe the upper part of the pipette with a soft-cloth or lint-free tissue soaked with distilled water or sterile water.

Single channel - Lower Part (Volumetric module)

The following components **only** can be immersed in a decontaminant solution: connecting nut, tip-ejector, tip-holder, piston assembly, seal and O-ring.

Multichannel - Lower Part (Volumetric module)

The following components **only** can be immersed in a cleaning solution: tip-ejector, ejector locks and ejector spacer.

- 1 Remove the tip-ejector and the ejector spacer.
- 2 Immerse the tip-ejector, ejector locks and ejector spacer in the cleaning solution or wipe them with a soft-cloth or lint-free tissue impregnated with the cleaning solution.
- 3 Rinse the components with distilled water.
- 4 Wipe the entire pipette with a soft cloth or lint-free tissue impregnated with the cleaning solution.
- 5 Wipe it with distilled water.
- 6 Leave the parts to dry by evaporation or wipe them with a clean soft-cloth or lint-free tissue.
- 7 Refit the tip-ejector as described in page 16 (Changing the tip-ejector).

13 - SPECIFICATIONS

PIPETMAN Neo is a high quality pipette that offers excellent accuracy and precision. The figures given in the "Gilson Maximum Permissible Errors" table (page 22) were obtained using PIPETMAN DIAMOND Tips. These figures are guaranteed only when genuine PIPETMAN DIAMOND Tips are used.

Each pipette is inspected and validated by qualified technicians in accordance with the Gilson Quality System. Gilson declares that its manufactured pipettes comply with the requirements of the ISO 8655 standard, by type testing.

The adjustment is carried out under strictly defined and monitored conditions (ISO 8655-6).

Gilson Maximum Permissible Errors

Single channel Model (Reference)	Volume (μL)		Maximum Permissible Errors			
			Gilson		ISO 8655	
			Systematic error (μL)	Random error (μL)	Systematic error (μL)	Random error (μL)
P2N (F144561) D10, DL10	Min.	0.2	± 0.024	≤ 0.012	± 0.08	≤ 0.04
		0.5	± 0.025	≤ 0.012	± 0.08	≤ 0.04
	Max.	2	± 0.030	≤ 0.014	± 0.08	≤ 0.04
P10N (F144562) D10, DL10	Min.	1	± 0.025	≤ 0.012	± 0.12	≤ 0.08
		5	± 0.075	≤ 0.030	± 0.12	≤ 0.08
	Max.	10	± 0.100	≤ 0.040	± 0.12	≤ 0.08
P20N (F144563) D200	Min.	2	± 0.10	≤ 0.03	± 0.20	≤ 0.10
		10	± 0.10	≤ 0.05	± 0.20	≤ 0.10
	Max.	20	± 0.20	≤ 0.06	± 0.20	≤ 0.10
P100N (F144564) D200	Min.	10	± 0.35	≤ 0.10	± 0.80	≤ 0.30
		50	± 0.40	≤ 0.12	± 0.80	≤ 0.30
	Max.	100	± 0.80	≤ 0.15	± 0.80	≤ 0.30
P200N (F144565) D200	Min.	20	± 0.50	≤ 0.20	± 1.60	≤ 0.60
		100	± 0.80	≤ 0.25	± 1.60	≤ 0.60
	Max.	200	± 1.60	≤ 0.30	± 1.60	≤ 0.60
P1000N (F144566) D1000	Min.	100	± 3	≤ 0.6	± 8	≤ 3
		500	± 4	≤ 1.0	± 8	≤ 3
	Max.	1000	± 8	≤ 1.5	± 8	≤ 3

Multichannel Model (Reference)	Volume (μL)		Maximum Permissible Errors			
			Gilson		ISO 8655	
			Systematic error (μL)	Random error (μL)	Systematic error (μL)	Random error (μL)
8x20 (F14401) and	Min.	2	± 0.10	≤ 0.08	± 0.40	≤ 0.20
		10	± 0.20	≤ 0.10	± 0.40	≤ 0.20
12x20 (F14402)	Max.	20	± 0.40	≤ 0.20	± 0.40	≤ 0.20
8x200 (F14403) and	Min.	20	± 0.50	≤ 0.25	± 3.2	≤ 1.2
		100	± 1.0	≤ 0.40	± 3.2	≤ 1.2
12x200 (F14404)	Max.	200	± 2.00	≤ 0.60	± 3.2	≤ 1.2



The data given in the tables conform to the ISO 8655-2 Standard.

With a precise pipetting technique (see Chapter 6 - General guidelines for good pipetting) the P2N model may be used to aspirate volumes as low as 0.1 μL and the P10N model as low as 0.5 μL .

14 - SPARE PARTS


Single channel models

Service Kit 1st level includes:


- 3 piston seals **C**
- 3 O-rings **D**
- 1 tip-holder **E**

Service Kit 2nd level includes:

- 1 push-button **A**
 - 1 connecting nut **B**
 - 1 tip-ejector **F**
- for the P2N and P10N models only:
- 1 tip-ejector **F1**
and 1 adapter **F2**




Multichannel models



F3

Tip-ejector (8 channels)




F4

Tip-ejector (12 channels)




G

Ejector lock



H1

Ejector spacer (8 channels)



H2

Ejector spacer (12 channels)

P2N (F144561) and P10N (F144562)

	Description	P2N	P10N
C+D+E	Service Kit 1st level	F144501	F144502
A+B+F1+F2	Service Kit 2nd level	F161970	F161971
C+D	Seal + O-ring (5 sets)	F144861	F144862
F2	Tip-ejector Adapter	F144879	F144879

P20N (F144563) and P100N (F144564)

	Description	P20N	P100N
C+D+E	Service Kit 1st level	F144495	F144496
A+B+F	Service Kit 2nd level	F161972	F161973
C+D	Seal + O-ring (5 sets)	F144863	F144864

P200N (F144565) and P1000N (F144566)

	Description	P200N	P1000N
C+D+E	Service Kit 1st level	F144497	F144498
A+B+F	Service Kit 2nd level	F161974	F161978
C+D	Seal + O-ring (5 sets)	F144865	F144866

8x20 (F14401) and 12x20 (F14402)

	Description	8x20	12x20
F3 - F4	Tip-ejector	F507005	F507006
G	Ejector lock	F507008	F507008
H1 - H2	Ejector spacer	F507001	F507003

8x200 (F14403) and 12x200 (F14404)

	Description	8x20	12x20
F3 - F4	Tip-ejector	F507005	F507006
G	Ejector lock	F507008	F507008
H1 - H2	Ejector spacer	F507001	F507003

EC DECLARATION OF CONFORMITY

The company,

GILSON S.A.S.

19, avenue des Entrepreneurs
BP. 145 - 95400 Villiers-le-Bel (France)
Tel.: +33(0) 1 34 29 50 00
Fax: +33(0) 1 34 29 50 20
www.gilson.com

Hereby certifies on its sole responsibility that the products listed below:

PIPETMAN Neo[®] Single channel
P2N, P10N, P20N, P100N, P200N, P1000N
PIPETMAN Neo[®] Multichannel
8x20, 12x20, 8x200, 12x200

comply with the requirements of the following European Directives:

98/79/EC*
on In Vitro Diagnostic Medical Devices

** Annex III, self-declared*

Villiers-le-Bel, July 2nd, 2012



A. El Sayed
General Manager



H. Ledorze
Quality Manager



WARRANTY

Gilson warrants this pipette against defects in material under normal use and service for a period of 12 months from the date of purchase.

This warranty shall not apply to pipettes which are subject to abnormal use and/or improper or inadequate maintenance (contrary to the recommendations given in the User's guide), including, but not limited to pipettes which have been subjected to physical damage, improper handling, or spillage or exposure to any corrosive environment. This warranty shall also be void in the event pipettes are altered or modified by any party other than Gilson or its designates. Gilson's sole liability under this warranty shall be limited to, at Gilson's sole option, repair or replacement of any defective components of pipettes or refund of the purchase price paid for such pipettes.

THE FOREGOING WARRANTY IS EXCLUSIVE AND GILSON HEREBY DISCLAIMS ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING ANY WARRANTIES OF MERCHANTABILITY AND ANY WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE, UNDER NO CIRCUMSTANCES SHALL GILSON BE LIABLE FOR ANY CONSEQUENTIAL, PUNITIVE, INDIRECT OR INCIDENTAL DAMAGES ARISING OUT OF ANY BREACH OF ANY EXPRESS OR IMPLIED WARRANTY.

LT801521/D - ©2012 Gilson SAS All rights reserved

July 2012 English
Printed in France

Specifications subject to change without notifications - errors omitted.

www.gilson.com

sales@gilson.com

service@gilson.com

training@gilson.com



Gilson, Inc. World Headquarters

3000 Parmenter Street | P.O. Box 620027 | Middleton, WI 53562-0027, USA
Tel: (1) 800-445-7661 or (1) 608-836-1551 | Fax: (1) 608-831-4451

Gilson S.A.S.

19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-Bel, France
Tel: (33-1) 34 29 50 00 | Fax: (33-1) 34 29 50 20